

## mRNA-EPO

mRNA-EPO is a ready-to-use mRNA product encoding human erythropoietin (EPO), designed for applications in hematology research, anemia model studies, mRNA drug development, and protein expression analysis.

Name	Cat.No	Scale	State
mRNA-EPO	HX-mRNA009	1mg	Lyophilized/Liquid

**Concentration:** 1mg/mL

**Storage Buffer:** ddH<sub>2</sub>O pH=7.0

**mRNA length:** 879 nt

**Amino acid sequence size:** 193 kDa

**Handling and Storage:** store at – 80°C for long term

**Amino acid sequence:**

MGVHECPAWLWLLLSLLSLPLGLPVLGAPPRLICDSRVLERYLLEAKEAENITGCAE  
HCSLNENITVPDTKVNIFYAWKRMEVGQQAVEVWQGLALLSEAVLRGQALLVNSSQP  
WEPLQLHVDKAVSGLRSLTLLRALGAQKEAISPPDAASAAPLRITADTFRKLFVYS  
NFLRGKCLKLYTGEACRTGDR